

# Anti-Phospho-AP2M1 (Thr156) Rabbit antibody



Catalog Number: 167473

## Product name

Anti-Phospho-AP2M1 (Thr156) Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, ICC/IF

## Antibody description

Phospho-AP2M1 (Thr156) Rabbit polyclonal antibody

## Dilutions

WB: 1/1000

ICC/IF: 1/20


## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr156 of human AP2M1

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of Phospho-AP2M1 (Thr156) in K562 cell lysates using Phospho-AP2M1 (Thr156) Rabbit pAb(1:1000 diluted). Predicted band size:50kDa. Observed band size:50kDa.

## Clonality

Polyclonal