

PI3K p110(beta) antibody



Catalog Number: 113873

Product name

PI3K p110(beta) antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human; other species not tested.

Validations

Antibody description

PI3K p110(beta) Rabbit Polyclonal antibody.
Positive WB detected in HepG2 cells, HeLa cells, Jurkat cells, K-562 cells, rat liver tissue. Positive IHC detected in human skin tissue, human lung cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 110 kDa




Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006219).
Purification method: Antigen affinity purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113873(PIK3CB antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human skin slide using Catalog No:113873(PIK3CB Antibody) at dilution of 1:50

Ig Type

Rabbit IgG



Applications

ELISA, WB, IF, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000


Immunohistochemistry of paraffin-embedded human skin slide using Catalog

PI3K p110(beta) antibody



Catalog Number: 113873

No:113873(PIK3CB Antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells, using PIK3CB antibody Catalog No:113873 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).