

RPL18A antibody

Catalog Number: 114873

Product name

RPL18A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RPL18A Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells. Positive IHC detected in human brain. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of RPL18A recombinant protein (Accession Number: NM_000980). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114873(RPL18A antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:114873(RPL18A antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:114873(RPL18A antibody) at dilution of 1:100 (under 40x lens)