

alpha Tubulin antibody



Catalog Number: 107873

Product name

alpha Tubulin antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

alpha Tubulin Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human testis tissue, human brain tissue, human kidney tissue, human medulloblastoma tissue, human ovary tissue, human placenta tissue, human skin tissue. Positive WB detected in HeLa cells, mouse cerebellum tissue. Observed molecular weight by Western-blot: 50-55 kDa

Validations



Preparation

This antibody was obtained by immunization of alpha Tubulin recombinant protein (Accession Number: NM_006009). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117300(tubulin-Alpha antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF, WB


Dilutions




Immunohistochemical of paraffin-embedded human testis using Catalog No:117300(tubulin-Alpha antibody) at dilution of 1:50 (under 10x lens)

alpha Tubulin antibody

Catalog Number: 107873



Immunohistochemical of paraffin-embedded human testis using Catalog No:117300(tubulin-Alpha antibody) at dilution of 1:50 (under 40x



lens)

Immunofluorescent analysis of HepG2 cells using Catalog No:117300(tubulin-Alpha Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG