

IRAK1 Polyclonal Antibody

Catalog Number: 161173

Product name

IRAK1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to IRAK1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 424-633 of human IRAK1 (NP_001020414.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:1000 - 1:3000


IF 1:50 - 1:200

Validations




Western blot - IRAK1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using IRAK1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunofluorescence - IRAK1 Polyclonal Antibody

Immunofluorescence analysis of C6 cells using IRAK1 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.


IRAK1 Polyclonal Antibody



Catalog Number: 161173

Immunofluorescence - IRAK1 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using IRAK1 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

Immunofluorescence - IRAK1 Polyclonal Antibody

Immunofluorescence analysis of NIH/3T3 cells using IRAK1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.