

CXCR1 Polyclonal Antibody



Catalog Number: 163173

Product name

CXCR1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CXCR1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 250-350 of human CXCR1 (NP_000625.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - CXCR1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CXCR1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.



Immunohistochemistry - CXCR1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using CXCR1 antibody at dilution of 1:100



(40x lens).


CXCR1 Polyclonal Antibody



Catalog Number: 163173

Immunohistochemistry - CXCR1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using CXCR1 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CXCR1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse spleen using CXCR1 antibody at dilution of 1:100 (40x lens).