

# NGFRAP1 antibody



Catalog Number: 113173

## Product name

NGFRAP1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NGFRAP1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 28kd

## Preparation

This antibody was obtained by immunization of NGFRAP1 recombinant protein (Accession Number: NM\_014380). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113173(NGFRAP1 antibody) at dilution of 1:400