

PRDM15 antibody

Catalog Number: 114173

Product name

PRDM15 antibody

Specificity

Human; other species not tested.

Antibody description

PRDM15 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells. Positive IP detected in MCF-7 cells. Observed molecular weight by Western-blot: 169 kDa

Preparation

This antibody was obtained by immunization of PRDM15 recombinant protein (Accession Number: XM_017028425). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114173(PRDM15 Antibody) at dilution of 1:300



IP Result of anti-PRDM15 (IP:Catalog No:114173, 5ug; Detection:Catalog No:114173 1:500) with MCF-7 cells lysate 2500ug.