

Product name

SGSH antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SGSH Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Transfected HEK-293 cells. Positive IHC detected in human pancreas cancer tissue, human liver cancer tissue.

Preparation

This antibody was obtained by immunization of SGSH recombinant protein (Accession Number: NM_000199). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115173(SGSH Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:115173(SGSH Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:115173(SGSH Antibody) at dilution of 1:200 (under 40x lens).