

# Arrestin C antibody

Catalog Number: 108173

## Product name

Arrestin C antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Arrestin C Rabbit Polyclonal antibody. Positive IHC detected in human retinoblastoma tissue. Positive WB detected in human brain tissue, mouse eye tissue, mouse heart tissue. Observed molecular weight by Western-blot: 43 kDa

## Preparation

This antibody was obtained by immunization of Arrestin C recombinant protein (Accession Number: BC012096). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108173 (ARR3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human retinoblastoma using Catalog No:108173 (ARR3 antibody) at dilution of 1:50 (under 10x lens)