

Anti-TMEM173 Rabbit antibody



Catalog Number: 168573

Product name

Rabbit IgG

Anti-TMEM173 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human

Dilutions

Antibody description

1:1000

TMEM173 Rabbit monoclonal antibody

1:20

Preparation

1:20

Antigen: Recombinant protein of human TMEM173

1:20

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of TMEM173 in K562,Rat Brain lysates using TMEM173 Rabbit mAb(1:1000 diluted).Predicted band size:42kDa.Observed band size:40kDa.