

# GIMAP4 antibody

Catalog Number: 110973

## Product name

GIMAP4 antibody

## Specificity

Human; other species not tested.

## Antibody description

GIMAP4 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 38 kDa

## Preparation

This antibody was obtained by immunization of GIMAP4 recombinant protein (Accession Number: NM\_018326). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110973(GIMAP4 antibody) at dilution of 1:1000



IP Result of anti-GIMAP4 (IP:Catalog No:110973, 3ug; Detection:Catalog No:110973 1:800) with Jurkat cells lysate 1080ug.