

Anti-Cystatin-B Rabbit antibody



Catalog Number: 167973

Product name

Anti-Cystatin-B Rabbit antibody

Specificity

Human, Rat

Antibody description

Cystatin-B Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Cystatin-B

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC


Dilutions

WB: 1/1000

IHC: 1/50

FC: 1/20

Validations



Western blot detection of Cystatin-B in HeLa,A549,HL-60 using Cystatin-B Rabbit mAb(1:1000 diluted)