

Human MDP1 (His tag) recombinant protein



Catalog Number: 518973

General Information

Gene Name Synonym

Magnesium-dependent phosphatase
1,MDP-1,EC=3.1.3.48,

Protein Construction

A DNA sequence encoding the Human MDP1
Met1-Ala176 is expressed with a 6His tag at the C-
terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM
TrisHCl,150mM NaCl,10%Glycerol,pH7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human MDP1 (His tag) recombinant protein