

HBQ1-Specific antibody



Catalog Number: 111273

Product name

HBQ1-Specific antibody

Specificity

Human; other species not tested.

Antibody description

HBQ1-Specific Rabbit Polyclonal antibody.
Positive WB detected in recombinant protein.
Observed molecular weight by Western-blot: 42kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005331).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



WB results of Catalog No:111273 (HBQ1 antibody) at dilution of 1:2000.