

# SLC5A6 Polyclonal Antibody



Catalog Number: 162273

## Product name

SLC5A6 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to SLC5A6

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human SLC5A6 (NP\_066918.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions


IHC 1:50 - 1:200

## Validations



Immunohistochemistry - SLC5A6 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver using SLC5A6 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - SLC5A6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using SLC5A6 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - SLC5A6 Polyclonal Antibody

# SLC5A6 Polyclonal Antibody



Catalog Number: 162273

---

Immunohistochemistry of paraffin-embedded

mouse heart using SLC5A6 antibody at dilution of 1:100 (40x lens).