

Anti-IL20RB antibody



Catalog Number: 100673

Product name

Anti-IL20RB antibody

Immunogen

[Rat IL20RB \(His Tag\) recombinant protein](#)

Specificity

Rat IL20RB

Antibody description

Mouse monoclonal to IL20RB

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Rat IL20RB (rR IL20RB; B8K1X9; Met1-Pro227). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS with 5% treha

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 10-20 μ g/ml

Validations



Lanes	A
Items	
Sample (whole cell lysate)	COLO205
Sample Volume (μ g/lane)	60
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	10-20 μ g/ml
Secondary Antibody	Dylight 800-labeled Antibody to Mouse IgG (H+L), at 1:2500 dilution.
	Developed using Odyssey imaging system.
Explanation	Predicted band size : 35 kDa Observed band size : 37 kDa Additional bands at : 85 kDa and 62 kDa (We are unsure as to the identity of these extra bands. They are possible non-specific bindings).

IL20RB Antibody, Mouse MAb, Western blot