

Rat VEGFR2/Flk-1/CD309/KDR (Fc Tag) recombinant protein



Catalog Number: 502673

General Information

Protein Construction

A DNA sequence encoding the rat KDR (AAH87029.1) (Met1-Glu760) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Activity

Immobilized mouse VEGF164 (Cat:502886) at 10 µg/ml (100 µl/well) can bind rat KDR-Fc, The EC_{50} of rat KDR-Fc is 0.07-0.15 µg/ml.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 20

Molecular Mass

The recombinant rat KDR comprises 982 amino acids and has a predicted molecular mass of 110 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat VEGFR2 / Flk-1 / CD309 / KDR Protein (Fc Tag) SDS-PAGE