

Mouse Growth Hormone Receptor / GHR / GHBP (His Tag) recombinant protein



Catalog Number: 503673

General Information

Gene Name Synonym

Somatotropin receptor; Growth hormone-binding protein; Serum-binding protein

Protein Construction

A DNA sequence encoding the extracellular domain of mouse GHR (NP_034414.2) (Met 1-Gln 273) was fused with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Thr 25

Molecular Mass

The recombinant mouse GHR consists of 260

amino acids and has a predicted molecular mass of 30.3 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rmGHR is approximately 40-45 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse Growth Hormone Receptor / GHR / GHBP Protein (His Tag) SDS-PAGE