

# REC8 antibody

Catalog Number: 114673

---

## Product name

REC8 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

REC8 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, K-562 cells. Observed molecular weight by Western-blot: 85-95kd

## Preparation

This antibody was obtained by immunization of REC8 recombinant protein (Accession Number: NM\_005132). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114673(REC8 antibody) at dilution of 1:1000