

# TST antibody

Catalog Number: 116373

## Product name

TST antibody

WB: 1:200-1:2000

IP: 1:200-1:1000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

TST Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in mouse kidney tissue, HEK-293 cells, HeLa cells, HepG2 cells, human kidney tissue, liver tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 33 kDa

## Validations




## Preparation

This antibody was obtained by immunization of TST recombinant protein (Accession Number: NM\_003312). Purification method: Antigen affinity purified.

WB result of anti-TST(Catalog No:116373) in Liver by Dr. Seow Chong Lee.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116373(TST antibody) at dilution of 1:50 (under 10x lens)

## Ig Type


Rabbit IgG

## Applications

ELISA, IHC, IP, WB

## Dilutions


Recommended Dilution:



# TST antibody

Catalog Number: 116373


Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116373(TST antibody) at dilution of 1:50



(under 40x lens)

mouse kidney tissue were subjected to SDS PAGE

followed by western blot with Catalog No:116373(TST Antibody) at dilution of 1:600



IP Result of anti-TST (IP:Catalog No:116373, 4ug; Detection:Catalog No:116373 1:300) with HEK-293 cells lysate 2800ug.