

# BBS6 antibody

Catalog Number: 108373

## Product name

BBS6 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BBS6 Rabbit Polyclonal antibody. Positive IF detected in hTERT-RPE1 cells. Positive WB detected in mouse testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 63kd

## Preparation

This antibody was obtained by immunization of BBS6 recombinant protein (Accession Number: NM\_018848). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

## Validations



IF result (cytoplasm and the base of cilia stain) of anti-BBS6 (Catalog No:108373; 1:50) with hTERT-RPE1 cell (MeOH fixed) by Dr. Moshe Kim.



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108373(BBS6 antibody) at dilution of 1:600