

MTFR1 antibody

Catalog Number: 112773

Product name

MTFR1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MTFR1 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue, A549 cells. Observed molecular weight by Western-blot: 37kd

Preparation

This antibody was obtained by immunization of MTFR1 recombinant protein (Accession Number: NM_014637). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112773(MTFR1 antibody) at dilution of 1:1000



IP Result of anti-MTFR1 (IP:Catalog No:112773, 4ug; Detection:Catalog No:112773 1:800) with mouse testis tissue lysate 4800ug.