

ROBO1-Specific antibody



Catalog Number: 114773

Product name

ROBO1-Specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ROBO1-Specific Rabbit Polyclonal antibody.
Positive WB detected in mouse brain tissue.
Observed molecular weight by Western-blot: 200 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002941).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114773(ROBO1-Specific antibody) at dilution of 1:1000