

H10N8 Neuraminidase (His Tag) recombinant protein



Catalog Number: 502083

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/duck/Guangdong/E1/2012(H10N8)) Neuraminidase (AFJ19054.1) (Gly38-Gly470) was expressed with a polyhistidine tag at the N-terminus.

Organism

H10N8

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant Neuraminidase of Influenza A virus (A/duck/Guangdong/E1/2012(H10N8))

consists 449 amino acids and predicts a molecular mass of 50.2 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H10N8 (A/duck/Guangdong/E1/2012) Neuraminidase / NA Protein (His Tag)