

# Anti-BRAP (1E7) Mouse antibody



Catalog Number: 168083

## Product name

Anti-BRAP (1E7) Mouse antibody

## Dilutions

WB: 1/1000;Other 1/500

## Specificity

Human

## Antibody description

BRAP (1E7) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human BRAP protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

WB

## Validations



Western blot detection of BRAP in A549 and HUVEC cell lysates using BRAP mouse mAb (1:500 diluted). Predicted band size: 67KDa. Observed band size: 67KDa.



Western blot analysis of extracts from CHO-K1 cells, transfected with a human pFLAG-CMV2-BRAP construct (A) or transfected with a human pFLAG-CMV2 construct (B), using BRAP mouse mAb (1:1000 diluted).