

Anti-UFD1L Rabbit antibody

Catalog Number: 169083

Product name

Anti-UFD1L Rabbit antibody

Specificity

Human, Mouse

Antibody description

UFD1L Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human UFD1L (NP_005650.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C . Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/500 - 1/2000


;IHC: 1/50 - 1/200

;IF: 1/50 - 1/100


Validations



Immunohistochemistry - UFD1L Polyclonal Antibody. Immunohistochemistry of paraffin-embedded human colon using UFD1L antibody at dilution of 1:100 (40x lens).



Immunofluorescence - UFD1L Polyclonal Antibody. Immunofluorescence analysis of MCF7



cells using UFD1L antibody .

Western blot - UFD1L Polyclonal Antibody. Western blot analysis of extracts of various cell lines, using UFD1L antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Exposure time:

Anti-UFD1L Rabbit antibody



Catalog Number: 169083

90s.