

Product name

HMGCL antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HMGCL Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human ovary tissue, human spleen tissue, human testis tissue. Positive WB detected in mouse liver tissue, human heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of HMGCL recombinant protein (Accession Number: NM_000191). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111483(HMGCL antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human liver using Catalog No:111483(HMGCL antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:111483(HMGCL antibody) at dilution of 1:100 (under 40x lens)