

Anti-CCL24 antibody



Catalog Number: 176483

Product name

Anti-CCL24 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CCL24

Preparation

This antigen of this antibody was synthetic peptide from human protein at aa range: 61-110

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

IHC-P

Dilutions

IHC-P: 1:50-200

Validations




Fig1:: Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200