

Mouse Tnfrsf9/Cd137/Ila/Ly63 (Fc tag) recombinant protein



Catalog Number: 519483

General Information

Gene Name Synonym

Secreted CD137 antigen □ Tumor necrosis factor receptor superfamily member 9 □ Tnfrsf9

Protein Construction

A DNA sequence encoding the Mouse Tnfrsf9/Cd137/Ila/Ly63 Val24-Leu187 is expressed with a Fc tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Tnfrsf9/Cd137/Ila/Ly63 (Fc tag)
recombinant protein