

GRIA3 Polyclonal Antibody



Catalog Number: 161883

Product name

GRIA3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to GRIA3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 200-490 of human GRIA3 (NP_000819.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:100

Validations



Immunofluorescence - GRIA3 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using GRIA3 antibody . Blue: DAPI for nuclear staining.