

Vopratelimab (Synonyms: ANTI-ICOS AGONIST MAB JTX-2011)



Catalog Number: 176883

Product name

Vopratelimab (Synonyms: ANTI-ICOS AGONIST MAB JTX-2011)

Specificity

Human ICOS / CD278

Antibody description

Anti-ICOS / CD278 Antibody (Vopratelimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SDS-PAGE

Anti-ICOS / CD278 Antibody (vopratelimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-ICOS / CD278 Antibody (vopratelimab) is 99.09%, determined by SEC-HPLC.