

# Anti-Cytokeratin 7 Rabbit antibody



Catalog Number: 168883

## Product name

Anti-Cytokeratin 7 Rabbit antibody

## Specificity

Human

## Antibody description

Cytokeratin 7 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Cytokeratin 7

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/2000-1/10000

ICC/IF: 1/50-1/200

FC: 1/20

IP: 1/20

## Validations



Western blot detection of Cytokeratin 7 in HepG2, MCF7 (negative for Keratin 7, as expected), HeLa cell lysates using Cytokeratin 7 Rabbit mAb (1:500 diluted). Predicted band size: 51 kDa. Observed band size: 51 kDa.