

# NSL1 Polyclonal Antibody



Catalog Number: 165183

---

## Product name

NSL1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to NSL1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-281 of human NSL1 (NP\_056286.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

## Validations



Immunofluorescence - NSL1 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using NSL1 antibody .