

# TMEM27 antibody



Catalog Number: 116183

## Product name

TMEM27 antibody

## Specificity

Human; other species not tested.

## Antibody description

TMEM27 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 23 kDa

## Preparation

This antibody was obtained by immunization of TMEM27 recombinant protein (Accession Number: NM\_020665). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116183(TMEM27 Antibody) at dilution of 1:300