

# CX3CL1 antibody



Catalog Number: 107183

## Product name

Monoclonal

CX3CL1 antibody

## Ig Type

## Specificity

Mouse IgG1

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

CX3CL1 Mouse Monoclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 90 kDa

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of CX3CL1 recombinant protein (Accession Number: NM\_002996). Purification method: Protein G purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



PC-3 cells were subjected to SDS PAGE followed by western blot with (CX3CL1 Antibody) at dilution of 1:1000