

CGB Polyclonal Antibody



Catalog Number: 163583

Product name

CGB Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CGB

Preparation

Antigen: Recombinant protein of human CGB

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - CGB Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CGB Antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.



Immunohistochemistry - CGB Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using CGB antibody at dilution of 1:100 (40x lens).