

# Anti-GLI3 antibody



Catalog Number: 176583

## Product name

Rabbit IgG

Anti-GLI3 antibody

## Applications

## Specificity

WB, ICC

Human

## Dilutions

## Antibody description

WB: 1:500-1:2,000

Rabbit monoclonal antibody to GLI3

ICC: 1:50-1:100

## Preparation

## Validations

This antigen of this antibody was recombinant protein within human gli3 aa 1261-1580 / 1580.

## Formulation

Liquid, 1\*TBS (pH7.4), 0.05% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type




Fig1.; Western blot analysis of Gli3 on PC-3M cell lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:200,000 dilution was used for 1 hour at room temperature.