

USHBP1 antibody

Catalog Number: 116583

Product name

USHBP1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse; other species not tested.

Antibody description

USHBP1 Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, human small intestine tissue. Positive IF detected in A549 cells. Positive WB detected in A549 cells, COLO 320 cells, mouse heart tissue. Observed molecular weight by Western-blot: 90-100 kDa

Preparation

This antibody was obtained by immunization of USHBP1 recombinant protein (Accession Number: NM_001321417). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116583(USHBP1 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:116583(USHBP1 Antibody) at dilution of 1:50 (under 40x lens)




Immunohistochemistry of paraffin-embedded

USHBP1 antibody

Catalog Number: 116583

human small intestine tissue slide using Catalog No:116583(USHBP1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (-20oc Ethanol) fixed A549 cells using Catalog No:116583(USHBP1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)