

# Anti-phospho-MEK4/MKK4 (Ser80) Rabbit antibody



Catalog Number: 167583

## Product name

Anti-phospho-MEK4/MKK4 (Ser80) Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

phospho-MEK4/MKK4 (Ser80) Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser80 of human MEK4/MKK4

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IP


## Dilutions

WB: 1/2000

IHC: 1/20

IP: 1/20

## Validations



Western blot detection of phospho-MEK4/MKK4 (Ser80) in A549, MCF-7, U251 cell lysates using phospho-MEK4/MKK4 (Ser80) Rabbit mAb(1:1000 diluted). Predicted band size: 44kDa. Observed band size: 44kDa.