

# MYO1B antibody

Catalog Number: 112983

---

## Product name

MYO1B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MYO1B Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, mouse brain tissue.

## Preparation

This antibody was obtained by immunization of MYO1B recombinant protein (Accession Number: NM\_001161819). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human heart using Catalog No:112983(MYO1B antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:112983(MYO1B antibody) at dilution of 1:50 (under 40x lens)