

TAF1C antibody



Catalog Number: 115983

Product name

TAF1C antibody

Specificity

human; other species not tested.

Antibody description

TAF1C Rabbit Polyclonal antibody. Positive WB detected in HT-1080 cells. Observed molecular weight by Western-blot: 95kd

Preparation

This antibody was obtained by immunization of TAF1C recombinant protein (Accession Number: NM_139353). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:115983(TAF1C antibody) at dilution of 1:300