

# Ziltivekimab (Synonyms: MEDI 5117 [WHO-DD])



Catalog Number: 176983

## Product name

Ziltivekimab (Synonyms: MEDI 5117 [WHO-DD])

## Specificity

Human IL-6 / IFN $\beta$ 2

## Antibody description

Anti-IL-6 / IFN $\beta$ 2 Antibody (Ziltivekimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

1 $\times$ PBS, pH 6.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-IL-6 / IFN $\beta$ 2 Antibody (ziltivekimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-IL-6 / IFN $\beta$ 2 Antibody (ziltivekimab) is more than 95%, determined by SEC-HPLC.