

# XBP1S antibody



Catalog Number: 116983

## Product name

Polyclonal

XBP1S antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

XBP1S Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells. Observed molecular weight by Western-blot: 55 kDa, 30 kDa

## Dilutions

Recommended Dilution:

## Preparation


WB: 1:200-1:1000

This antibody was obtained by immunization of XBP1S recombinant protein (Accession Number: NM\_005080). Purification method: Antigen affinity purified.

## Validations

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:116983(XBP1 Antibody) at dilution of 1:300

## Clonality