

# XBP1S antibody



Catalog Number: 116983

## Product name

XBP1S antibody

## Specificity

Human; other species not tested.

## Antibody description

XBP1S Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells. Observed molecular weight by Western-blot: 55 kDa, 30 kDa

## Preparation

This antibody was obtained by immunization of XBP1S recombinant protein (Accession Number: NM\_005080). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:116983(XBP1 Antibody) at dilution of 1:300