

ALK,CD246 antibody

Catalog Number: 107983

Product name

ALK,CD246 antibody

Immunogen

[Human ALK,CD246 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ALK,CD246 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue, human brain tissue.

Preparation

This antibody was obtained by immunization of ALK,CD246 recombinant protein (Accession Number: NM_004304). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:
IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:107983(ALK,CD246 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:107983(ALK,CD246 Antibody) at dilution of 1:50 (under 40x lens)