

LMO2 antibody



Catalog Number: 112283

Product name

LMO2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

LMO2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 17-25 kDa

Preparation

This antibody was obtained by immunization of LMO2 recombinant protein (Accession Number: XM_017017731). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112283(LMO2 Antibody) at dilution of 1:300