

# Anti-RAP2C antibody

Catalog Number: 176283

---

**Product name**

Polyclonal

Anti-RAP2C antibody

**Ig Type****Specificity**

Rabbit IgG

Human, Mouse

**Applications****Antibody description**

WB

Rabbit polyclonal antibody to RAP2C

**Dilutions****Preparation**

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rap2c. aa range:97-146

**Validations****Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


**Clonality**

Fig1.; Western Blot analysis of various cells using Rap 2C Polyclonal Antibody diluted at 1: 500