

Canine CTLA-4/CD152/CD152 (Fc Tag) recombinant protein



Catalog Number: 503683

General Information

Gene Name Synonym

Cytotoxic T-lymphocyte protein 4

Protein Construction

A DNA sequence encoding the canine CD152 (NP_001003106.1) (Met1-Asp161) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Canine

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Lys 36

Molecular Mass

The recombinant canine CD152 consists 364 amino acids and predicts a molecular mass of 40.1 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



CD152 protein SDS-PAGE