

STARS/ABRA antibody



Catalog Number: 115683

Product name

STARS/ABRA antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

STARS/ABRA Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human heart tissue. Positive WB detected in mouse heart tissue, human heart tissue. Observed molecular weight by Western-blot: 45 kDa




Preparation

This antibody was obtained by immunization of STARS/ABRA recombinant protein (Accession Number: NM_139166). Purification method: Antigen affinity purified.

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115683(ABRA antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human heart using Catalog No:115683(ABRA antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IF, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human heart using Catalog No:115683(ABRA

STARS/ABRA antibody

Catalog Number: 115683

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:115683(ABRA Antibody) at dilution of 1:50 and Rhodamine-labeled goat anti-rabbit IgG (red).